



# Recombinant D-tagatose-1,6-bisphosphate aldolase subunit kbaZ (kbaZ)

<b>Product Code</b>	CSB-EP812540EGX
<b>Abbreviation</b>	kbaZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FDB3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKHLTEMVRQ HKAGKTNGIY AVCSAHPLVL EAAIRYASAN QTPLLIATS NQVDQFGGYT GMTPADFRGF VCQLADSLNF PQDALILGGD HLGPNRWQNL PAAQAMANAD DLIKSVAAG FKKIHLDCSM SCQDDPIPLT DDIVAERAAR LAKVAEETCL EHFGEADLEY VIGTEVPVPG GAHETLSELA VTTPDAARAT LEAHRHAFEK QGLNAIWPRI IALVVQPGVE FDHTNVIDYQ PAKATXLSQM XENYETLIFE AHSTDYQTPQ SLRQLVIDHF AILKVGPAIT FALREALFSL AAIEEELVPA KACSGLRQVL ENVMLDRPEY WQSHYHGDGN ARRLARGYSY SDRVRYYPWD SQIDDAFAHL VRNLADSPIP LPLISQYLPL QYVKVRSSEL QPTPRELIIN HIQDILAQYH TACEGQ
<b>Source</b>	E.coli
<b>Target Names</b>	kbaZ
<b>Protein Names</b>	Recommended name: D-tagatose-1,6-bisphosphate aldolase subunit kbaZ
<b>Expression Region</b>	1-426
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.