



# Recombinant Chaperone protein hscA (hscA)

<b>Product Code</b>	CSB-BP812548EGX
<b>Abbreviation</b>	hscA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FF44
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALLQISEPG LSAAPHQRRRL AAGIDLGTTN SLVATVRSQG AETLADHEGR HLLPSVVHYQ QQGHSVGYDA RTNAALDTAN TISSVKRLMG RSLADIQQRY PHLPYQFQAS ENGLPMIETA AGLLNPRVRS ADILKALAAR ATEALAGELD GVVITVPAYF DDAQRQGTKD AARLAGLHVL RLLNEPTAAA IAYGLDSGQE GVIADVYDLGG GTFDISILRL SRGVFEVLAT GGDSALGGDD FDHLLADYIR EQAGIPDRSD NRVQRELLDA AIAAKIALSD ADSVTVNVAG WQGEISREQF NELIAPLVKR TLLACRRALK DAGVEADEVL EVVVVGGSTR VPLVRERVGE FFGRPPLTSI DPKVVAIGA AIQADILVGN KPDSEMLLLD VIPLSLGLET MGGLVEKVIP RNTTIPVARA QDFTTFKDGQ TAMSIHVMQG ERELVDQCRS LARFALRGIP ALPAGGAHIR VTFQVDADGL LSVTAMEKST GVEASIQVKP SYGLTDSEIA SMIKDSMSYA EQDVKARMLA EQKVEAARVL ESLHGALAAD AALLSAAERQ AIDDAAAHLS EVAQGDDVDA IEQAIKNVDK QTQDFAARRM DQSVRRALKG HSVDEV
<b>Source</b>	Baculovirus
<b>Target Names</b>	hscA
<b>Protein Names</b>	Recommended name: Chaperone protein hscA Alternative name(s): Hsc66
<b>Expression Region</b>	1-616
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.