



Recombinant Chaperone protein hscA (hscA)

Product Code	CSB-EP812548EGX
Abbreviation	hscA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FF44
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MALLQISEPG LSAAPHQRRRL AAGIDLGTTN SLVATVRSGQ AETLADHEGR HLLPSVVHYQ QQGHSVGYDA RTNAALDTAN TISSVKRLMG RSLADIQQRY PHLPYQFQAS ENGLPMIETA AGLLNPRVRS ADILKALAAR ATEALAGELD GVVITVPAYF DDAQRQGTKD AARLAGLHVL RLLNEPTAAA IAYGLDSGQE GVIAVYDLGG GTFDISILRL SRGVFEVLAT GGDSALGGDD FDHLLADYIR EQAGIPDRSD NRVQRELLDA AIAAKIALSD ADSVTVNVAG WQGEISREQF NELIAPLVKR TLLACRRALK DAGVEADEVL EVVVVGGSTR VPLVRERVGE FFGRPPLTSI DPKVVAIGA AIQADILVGN KPDSEMLLLD VIPLSLGLET MGGLVEKVIP RNTTIPVARA QDFTTFKDGQ TAMSIHVMQG ERELVDQCRS LARFALRGIP ALPAGGAHIR VTFQVDADGL LSVTAMEKST GVEASIQVKP SYGLTDSEIA SMIKDSMSYA EQDVKARMLA EQKVEAARVL ESLHGALAAD AALLSAAERQ AIDDAAAHLS EVAQGDDVDA IEQAIKNVDK QTQDFAARRM DQSVRRALKG HSVDEV
Source	E.coli
Target Names	hscA
Protein Names	Recommended name: Chaperone protein hscA Alternative name(s): Hsc66
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.