



Recombinant Anaerobic glycerol-3-phosphate dehydrogenase subunit B (glpB)

Product Code	CSB-EP812553EGX-B
Abbreviation	glpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FFN4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRFDTVIMGG GLAGLLCGLQ LQKHGLRCAI VTRGQSALHF SSGSLDLLSH LPDGQPVTDI HSGLESLRQQ APAHPYTLLG PQRVLDLACQ AQALIAESGA QLQGSVELAH QRITPLGTLR STWLSSPEVP VWPLPAKKIC VVGISGLMDF QAHLAAASLR ELDLAVETA E IELPELDVLR NNATEFRAVN IARFLDNEEN WPLIIDALIP VANTCEMILM PACFGLADDK LWRWLNEKLP CSLMLLPTLP PSVLGIRLQN QLQRQFVRQG GVWMPGDEVK KVTCKNGVVN EIWTRNHADI PLRPRFVULA SGSFFSGGLV AERDGI REPI LGLDVLQTAT RGEWYKGDFF APQPWQQFGV TTDEALRPSQ AGQT IENLFA IGSVLGGFDP IAQGCGGGVC AVSALHVAQQ IAQRAGGQQ
Source	E.coli
Target Names	glpB
Protein Names	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit B Short name= Anaerobic G-3-P dehydrogenase subunit B Short name= Anaerobic G3Pdhase B EC= 1.1.5.3
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.