



Recombinant Arabidopsis thaliana 7-hydroxymethyl chlorophyll a reductase, chloroplastic (HCAR)

| | |
|--------------------------|---|
| Product Code | CSB-YP812734DOA |
| Abbreviation | HCAR |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8GS60 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SSSSRSKDMN LEPKKKVKLR EDWREKSRPI PPGGTYPAKD HCSQCGLCDT YYIAHVKEAC AFLGDGMSRI ESLEPVVHGR GRKADSLQDT YFGVHQEQLY ARKLPVEGA QWTGIVTTIA IEMLKSNMVE AVVCVQSDPE DRLSPRPVLA RTPEEVLAAR GVKPTLSPNL NTLELIEASG VKRLLFCGVG CQVQALRSVE QHLNLEKLYV LGTNCVDNGT RDGLDKFLKA ASKEPETVLH YEFMQDYKVV LKHLDGHIIE VPFYSLPAND LVDVIAPSCY SCFDYTNALA DLVIGYMGVP KYSGLNMTDH PQYITVRNER GKEMLSLVEN LLEITPTISS GDRRPFVTET VKADDAAKFG QGPAQPAPLF VGNIIAFILN LVGPKGLEFA RYSLDYHTIR NYLYVNRKWG KQRANTHMPS YAKKIVEMYN KNGQIDKMLS KK |
| Source | Yeast |
| Target Names | HCAR |
| Protein Names | Recommended name: 7-hydroxymethyl chlorophyll a reductase, chloroplastic EC= 1.17.7.2 |
| Expression Region | 21-462 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.