



Recombinant Arabidopsis thaliana B3 domain-containing protein REM7 (REM7)

Product Code	CSB-MP812772DOA
Abbreviation	REM7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H2D0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MADPVLQSPK NSHFFQPLLP GFDSYLNIPV KFFSKRIQGR NEGRTVELRT DASEKTWQVK IQGRRLTVGW KEFASAHDRF VGDIIVFRHE GSLVFHVTAL GPSCCEIQYV PSCNNQENIS NLSMKQSIKT ELESSLDEDK VNMGKFPRKK HVKKRIPEAE AKSFSSDKSC FVAHVTDNSL REDTLFLPRK FDRSDGLIKG SNKIVLMNEE ARTWTLILKF RNSSRTFYMR GGWTSFMSM GLNLEIPSRL N
Source	Mammalian cell
Target Names	REM7
Protein Names	Recommended name: B3 domain-containing protein REM7 Alternative name(s): Auxin response factor 25 Protein REPRODUCTIVE MERISTEM 7
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.