



Recombinant *Oryza sativa* subsp. japonica Tubby-like F-box protein 11 (TULP11)

Product Code	CSB-BP812774OFG
Abbreviation	TULP11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H485
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSFRSVIQEV KGEIGAISRR GFRSRPGRVR RVAAAAEPP DESSAAALVM RESCWTQLPP ELLREVLARV EESEGWWPRR RDVVACAGVC RSWRGIVREI VRTPEASGNL TFPISLKQPG PRDAPMKCFI VRNRTTQTY Y LYIGLTDALT DDGKFLLAAR KCRRTTCTEY LISLDMNDIS KRTDSYVGKL RSNFLGTKFT IYDAHPPYAG DVISKGQSAR VIGSNHLSR IPAGNYPVSH ISYELNVLGS RGPRRMHCAM DSIPVSAIEQ GGTAPTQTEF PLSYHESFTS IPFFKSKSVR ANNSTASLLT QNGSKLV LKN KSPRWHEHLQ CWCLNFHGRV TVASVKNFQL VASDESNPTN QEHDDVILQF GKVGKDMFTM DYRYPISAFQ AFAICLSSFD TKIACE
Source	Baculovirus
Target Names	TULP11
Protein Names	Recommended name: Tubby-like F-box protein 11 Short name= OsTLP11 Alternative name(s): Tubby-like F-box protein 7 Short name= OsTLP7
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.