



Recombinant Human Tumor protein p53-inducible nuclear protein 2 (TP53INP2)

Product Code	CSB-EP812888HU-B
Abbreviation	TP53INP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXH6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFQRLSSLFF STPSPPEDPD CPRAFVSEED EVDGWLII DL PDSYAAPPSP GAAPAPAGRP PPAPSLMDES WFVTPPACFT AEGPGLGPAR LQSSPLEDLL IEHPMSVYV TGSTIVLEPG SPSPLPDAAL PDGDLSEGEL TPARREPRAA RHAAPL PARA ALLEKAGQVR RLQRRARQRAE RHALSAKAVQ RQNRARES RP RRSKNQSSFI YQPCQRQFN Y
Source	E.coli
Target Names	TP53INP2
Protein Names	Recommended name: Tumor protein p53-inducible nuclear protein 2 Alternative name(s): p53-inducible protein U Short name= PIG-U
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.