



Recombinant Human Zinc finger protein DZIP1L (DZIP1L)

Product Code	CSB-MP812905HU
Abbreviation	DZIP1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYY4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQSPAATAEG LSGPLFGAYT FPTFKFQPRH DSMDWRRIST LDVDRVAREL DVATLQENIA GITFCNLDRE VCSRCGQPVD PALLKVLRLA QLIIEYLLHC QDCLSASVAQ LEARLQTSLG QQQRGQQELG RQADELKGVR EESRRRRKMI STLQQLLMQT GTHSYHTCHL CDKTFMNATF LRGHIQRRHA GVAEGGKQKK QEQPVVEEVLE ELRAKLKWTQ GELEAQREAE RQRQLQEAE IHQREIEAKK EFDKWKEQEW TKLYGEIDKL KKLFWDEFKN VAKQNSTLEE KLRALQSHSV MESKLGSLRD EESEEWLRQA RELQALREKT EIQKTEWKRK VKELHEEHMA EKKELQEENQ RLQASLSQDQ KKAQAQSQSQ ISTLRAQLQE QARIIASQEE MIQSLSLRKV EGIHKVPAKAV DTEEDSPEEE MEDSQDEQHK VLAALRRNPT LLKHFRPILE DTLEEKLESM GIRKDAKGIS IQTLRHLESL LRVQREQKAR KFSEFLSLRG KLVKEVTSRA KERQENGAVV SQPDGQPSVK SQQSTLVTRE AQPKTRTLQV ALPSTPAEPP PPTRQSHGSH GSSLTQVSAP APRPGLHGPS STPPSSGPGM STPPFSSEED SEGDRVQRVS LQPPKVPSRM VPRPKDDWDW SDTETSEENA QPPGQGSGTL VQSMVKNLEK QLEAPAKKPA GGVSLFFMPN AGPQRAATPG RKPQLSEDES DLEISSLEDL PLDLQREKP KPLSRSKLPE KFGTGPQSSG QPRVPAW
Source	Mammalian cell
Target Names	DZIP1L
Protein Names	Recommended name: Zinc finger protein DZIP1L Alternative name(s): DAZ-interacting protein 1-like protein
Expression Region	1-767
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.