



# Recombinant Mouse WAS/WASL-interacting protein family member 1 (Wipf1)

<b>Product Code</b>	CSB-MP812968MO
<b>Abbreviation</b>	Wipf1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K117
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPVPPPPAPP PPPTFALANT EKPTLNKTEQ AGRNALLSDI SKGKKLKKTV TNDRSAPILD KPKGAGASAG GYGGGGGGGG GGGGGGGGSG GNFGGGGPPG LGGLFQAGMP KLRSTANRDN DSGGSRPPIL PPGGRATSAK PFSPPSGPGR FPAPSPGHRG GPPEPPRNRM PPRPDVGSK PDSLPPVVPN TPRVPSSLH NRGSPAGLGA PRPPFPGNRG AAFGAGSARQ NPSGSSSPFP RPPLPPTPSR ALDDKPPPPP PPVGNRPSMH REAVPPPPSQ TSKPPVPSTP RPGLGSQAPP PPPPSRPGP PPLPPASNDE IPRLPQRNLS LTSSAPPLPS PGRSGPLPPP PSERPPPPVR DPPGRSGPLP PPPPINRNGS TARALPATPQ LPSRSGMDSP RSGPRPPLPP DRPGAGAPPP PPPSTSVRNG FQDSSCEDEW ESRFYFHPIS DLPPPEPYVP TTKTYPKLA RNESRSGSNR RERGAPPLPP IPR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Wipf1
<b>Protein Names</b>	Recommended name: WAS/WASL-interacting protein family member 1 Alternative name(s): Wiskott-Aldrich syndrome protein-interacting protein Short name= WASP-interacting protein
<b>Expression Region</b>	1-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a protein that plays an important role in the organization of the actin cytoskeleton. The encoded protein binds to a region of Wiskott-Aldrich syndrome protein that is frequently mutated in Wiskott-Aldrich syndrome, an X-linked recessive disorder. Impairment of the interaction between these two proteins may contribute to the disease. Two transcript variants encoding the same protein have been identified for this gene.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.