



Recombinant Non-canonical purine NTP phosphatase (yjjX)

Product Code	CSB-EP813840EGX
Abbreviation	yjjX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FA41
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MHQVVCATTN PAKIQAILQA FHEIFGEGSC HIASVAVESG VPEQPFGSEE TRAGARNRVA NARRLLPEAD FWVAIEAGID GDSTFSWVVI ENTSQRGEAR SATLPLPAVI LEKVREGEAL GPVMSRYTGI DEIGRKEGAI GVFTAGKLTR TSVYHQAVIL ALSPFHNAVY QPLQA
Source	E.coli
Target Names	yjjX
Protein Names	Recommended name: Non-canonical purine NTP phosphatase EC= 3.6.1.- Alternative name(s): Non-standard purine NTP phosphatase Nucleoside-triphosphate phosphatase Short name= NTPase
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.