



Recombinant Ornithine carbamoyltransferase, catabolic (arcB)

Product Code	CSB-MP813842EGX
Abbreviation	arcB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FAE1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MATSLKNRNF LKLLDYTPAE IQYLIDLAIN LKAAKKSNE KQTLVGKNIA LIFEKSSTRT RCAFEEAAFD QGAQVTYIGP SGSQIGHKES MKDTARVLGR MYDGEYRGY GQNIVEELGE FAGVPVWNGL TNEFHPTQIL ADLMTMLEHA PGKTLPELSF AYLGDARNNM GNSLMVGA AKMGMDIRLVAP KSFWPDEALV TQCREIASVT GARITLTEDV EEGVYDVDFL YTDVWVSMGE PKEAWAERVS LMAPYQINQQ VITATRNPEV KFMHCLPAFH NEHTTVGREI EMAYGLKGLE VTEEVFESAH SIVFDEAENR MHTIKAVMVA TLGD
Source	Mammalian cell
Target Names	arcB
Protein Names	Recommended name: Ornithine carbamoyltransferase, catabolic Short name= OTCase EC= 2.1.3.3
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.