



Recombinant L-ribulose-5-phosphate 3-epimerase UlaE (ulaE)

Product Code	CSB-MP813843EGX
Abbreviation	ulaE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FAI7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLSKQIPLGI YEKALPAGEC WLERLQLAKT LGDFDFVEMSV DETDDRLSRL DWSREQRLAL VNAIVETGVR VPSMCLSAHR RFPLGSEDDA VRAQGLEIMR KAIQFAQDVG IRVIQLAGYD VYYQEANNET RRRFRDGLKE SVEMASRAQV TLAMEIMDYP LMNSISKALG YAHYLNPNWF QLYPDIGNLS AWDNDVQMEL QAGIGHIVAV HVKDTKPGVF KNVPGEGVV DFERCFETLK QSGYCGPYLI EMWSETAEDP AAEVAKARDW VKARMAKAGM VEEA
Source	Mammalian cell
Target Names	ulaE
Protein Names	Recommended name: L-ribulose-5-phosphate 3-epimerase UlaE EC= 5.1.3.22 Alternative name(s): L-ascorbate utilization protein E L-xylulose-5-phosphate 3-epimerase
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.