



Recombinant Molybdenum cofactor guanylyltransferase (mobA)

Product Code	CSB-MP813850EGX
Abbreviation	mobA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBH7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNLMTMITGV VLAGGKARRM GGVDKGLLEL NGKPLWQYVA DALMTQLSHV VNNANRHQEI YQVSGLKVIE DSLADYPGPL AGMLSVMQQE AGEWFLFCPC DTPYIPHDLA ARLTHQRKDA PVVWVHDGER DHPTIALVNR AIEPLLEYL QAGERRVMAF MRLAGGHAVD FSDRKEAFIN VNTPEELARW QEKR
Source	Mammalian cell
Target Names	mobA
Protein Names	Recommended name: Molybdenum cofactor guanylyltransferase Short name= MoCo guanylyltransferase EC= 2.7.7.77 Alternative name(s): GTP:molybdopterin guanylyltransferase Mo-MPT guanylyltransferase Molybdopterin guanylyltransferase
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.