



Recombinant Guanosine-5'-triphosphate,3'-diphosphate pyrophosphatase (gppA)

Product Code	CSB-MP813853EGX
Abbreviation	gppA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBR0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MGSTSSLYAA IDLGSNSFHM LVVREVAGSI QTLTRIKRKV RLAAGLNSEN ALSNEAMERG WQCLRLFAER LQDIPPSQIR VVATATLRLA VNAGDFIAKA QEILGCPVQV ISGEEEARLI YQGVAHTTGG ADQRLVVDIG GASTELVTGT GAQTTSLFSL SMGCVTWLER YFADRNLGQE NFDAAEKAAR EVLRPVADEL RYHGWKVCVG ASGTVQALQE IMMAQGMDER ITLEKLQQLK QRAIHCGRLE ELEIDGLTLE RALVFPSTGLA ILIAIFTELN IQCMTLAGGA LREGLVYGML HLTVEQDIRS RTLRNIQRRF MIDIDQAQRV AKVAANFYDQ VENEWHLEAI SRDLLISACQ LHEIGLSVDF KQAPQHAAYL VRNLDLPGFT PAQKLLATL LLNQTNPVDL SSLHQQNAV PPRVAEQLCRL LRLAIIFASR RRDDLVPENT LQANHELLTL TLPQGWLQHQH PLGKEIIDQE SQWQSYVHWP LEVH
Source	Mammalian cell
Target Names	gppA
Protein Names	Recommended name: Guanosine-5'-triphosphate,3'-diphosphate pyrophosphatase EC= 3.6.1.40 Alternative name(s): Guanosine pentaphosphate phosphohydrolase pppGpp-5'-phosphohydrolase
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.