



Recombinant HTH-type transcriptional regulator yidZ (yidZ)

Product Code	CSB-MP813855EGX
Abbreviation	yidZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBU9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKKSITLDDL NLLLCLQLLM QERSVTKAAK RMNVTPSAVS KSLAKLRAWF DDPLFVNSPL GLSPTPLMVS MEQNLAEWMQ MSNLLLDKPH HQTPRGLKFE LAAESPLMMI MLNALSRIY QRYPQATIKL RNWDYDSLDA ITRGEVDIGF SGRESHPRSR ELLSSLPLAI DYEVLFSVDP CVWLRKDHPA LHEAWNLDTF LRYPHISICW EQSDTWALDN VLQELGRERT IAMSLPEFEQ SLFMAAQPDN LLLATAPRYC QYYNQLHQLP LVALPLPFDE SQQKKLEVPF TLLWHKRNSR NPKIVWLRET IKNLYASMA
Source	Mammalian cell
Target Names	yidZ
Protein Names	Recommended name: HTH-type transcriptional regulator yidZ
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.