



# Recombinant Siroheme synthase (cysG)

<b>Product Code</b>	CSB-MP813862EGX
<b>Abbreviation</b>	cysG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FCW8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDHLPIFCQL RDRDCLIVGG GDVAERKARL LLDAGARLTV NALAFIPQFT AWADAGMLTL VEGPFDESLL DTCWLAI AAT DDDALNQRVS EAAEARRIFC NVVDAPKAAS FIMPSIIDRS PLMVAVSSGG TSPVLARLLR EKLESLLPLH LGQVAKYAGQ LRGRVKQQFA TMGERRRFWE KLFVNDRLAQ SLANNDQKAI TETTEQLINE PLDHRGEVVL VGAGPGDAGL LTLKGLQQIQ QADV VVYDRL VSDDIMNLVR RDADRVFVGK RAGYHCVPQE EINQILLREA QKGKRVVRLK GGDPFIFGRG GEELETLCNA GIPFSVVPGI TAASGCSAYS GIPLTHR DYA QSVRLITGHL KTG GELD WEN LAAEKQTLVF YMGLNQAATI QK K LIEY GMP GEMPVAIVEN GTAVTQRVID GTLTQLGELA QQMNSPSLII IGRVVGLRDK LNWFSNH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	cysG
<b>Protein Names</b>	Recommended name: Siroheme synthase Including the following 3 domains: Uroporphyrinogen-III C-methyltransferase Short name= Urogen III methylase EC= 2.1.1.107 Alternative name(s): SUMT Uroporphyrinogen III methylase S
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.