



# Recombinant UPF0597 protein yhaM (yhaM)

<b>Product Code</b>	CSB-MP813867EGX
<b>Abbreviation</b>	yhaM
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FDC6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFDSTLNPLW QRYILAVQEE VKPALGCTEP ISLALAAVA AAELEGPVER VEAWVSPNLM KNGLGVTVPG TGMVGLPIAA ALGALGGNAN AGLEVLKDAT AQAIADAKAL LAAGKVSVKI QEPCDEILFS RAKVWNGEKW ACVTIVGGHT NIVHIETHNG VVFTQQACVT EGEQESPLTV LSRTTLAEIL KFVNEVPFAA IRFILDSAKL NCALSQEGLS GNWGLHIGAT LEKQCARGLL AKDLSSSIVI RTSAASDARM GGATLPAMSN SGSGNQGITA TMPVVVVAEH FGADDERLAR ALMLSHLSAI YIHNQLPRLS ALCAATTAAM GAAAGMAWLV DGRYETISMA ISSMIGDVSG MICDGASNSC AMKVSTSASA AWKAVLMALD DTAVTGNEGI VAHDVEQSIA NLCALASHSM QQTDRQIIEI MASKAR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	yhaM
<b>Protein Names</b>	Recommended name: UPF0597 protein yhaM
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.