



Recombinant Arabinose 5-phosphate isomerase KpsF (kpsF)

Product Code	CSB-BP813870EGX
Abbreviation	kpsF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FDQ2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSERHLPDDQ SSTIDPYLIT SVRQTLAEQS AALQNLSKQL DSGQYQRVLN LIMNCKGHVI LSGMGKSGHV GRKISATLAS TGTPSFFIHP AEA FHGDLGM ITPYDLLILI SASGETDEIL KLVPSLKNFG NRIIAITNNG NSTLAKNADA VLELHMANET CPNNLAPPTS TTLTMAIGDA LAIAMIHQRK FMPNDFARYH PGGSLGRLL TRVADVMQHD VPAVQLDASF KTVIQRITSG CQGMVMVEDA EGGLAGIITD GDLRRFMEKE GSLTSATAAQ MMTREPLTLP EDTMIIEAEE KMQKHRVSTL LVTNKANKVT GLVRIFD
Source	Baculovirus
Target Names	kpsF
Protein Names	Recommended name: Arabinose 5-phosphate isomerase KpsF Short name= API EC= 5.3.1.13 Alternative name(s): K-antigen-specific arabinose 5- phosphate isomerase Short name= K-API
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.