



# Recombinant Ribosomal RNA small subunit methyltransferase F (rsmF)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP813893EGX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Abbreviation</b>      | rsmF                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                      |
| <b>Uniprot No.</b>       | Q8FGS9                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Sequence</b>          | MAQHTVYFPD AFLTQMREAM PSTLSFDDFL AACQRPLRRS IRVNTLKTSV<br>ADFLQLTAPY GWTLTPIPWC EEGFWIERDN EDALPLGSTA EHLSGLFYIQ<br>EASSMLPVAA LFADGNAPQR VMDVAAAPGS KTTQIAARMN NEGAILANEF<br>SASRVKVLHA NISRCGISNV ALTHFDGRVF GAAVPEMFDA ILLDAPCSGE<br>GVVRKDPDAL KNWSPESNQE IAATQRELID SAFHALRPGG TLVYSTCTLN<br>REENEAVCLW LKETYHDAVE FLPLGDLFPG ANKALTEDGF LHVFPQIYDC<br>EGFFVARLRK TQAIPVLPAP KYKVGNFPS PVKDREAGQI RQAAASVGLN<br>WDENLRLWQR DKELWLFPVG IEALIGKVRF SRLGIKLAET HNKGYRWQHE<br>AVIALASPDN VNAFELTPQE AEEWYRGRDV YPQAAPVADD VLVTFQHQP<br>GLAKRIGSRL KNSYPRELVR DGKLFTGNA |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Names</b>      | rsmF                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Names</b>     | Recommended name: Ribosomal RNA small subunit methyltransferase F EC= 2.1.1.178 Alternative name(s): 16S rRNA m5C1407 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RsmF                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-479                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.                                                                                                                                                                                                                                                                                                                                                                                                                              |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.