



Recombinant Glutamine--tRNA ligase (glnS)

| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP813913EGX |
| Abbreviation | glnS |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8FJW4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | SEAEARPTN FIRQIIDDL ASGKHTTVHT RFPPEPNGYL HIGHAKSICL NFGIAQDYKG QCNLRFD DDTN PVKEDIEYVD SIKNDVEWLG FHWSGNVRYS SDYFDQLHAY AIELINKGLA YVDELTPEQI REYRGTLTQP GKNSPYRDRS VEENLALFEK MRTGGFEEGK ACLRAKIDMA SPFIVMRDPV LYRIKFAEHH QTGNKWCYIP MYDFTHCISD ALEGITHSLC TLEFQDNRRLL YDWVLDNITI PVHPRQYEFS RLNLEYTVMS KRKLNLLVTD KHVEGWDDPR MPTISGLRRR GYTAASIREF CKRIGVTKQD NTIEMASLES CIREDLNENA PRAMAVIDPV KLVIENYQGE GEMVTMPNHP NKPEMGSRQV PFSGEIWIDR ADFREEANKQ YKRLVLGKEV RLRNAYVIKA ERVEKDAEGN ITTIFCTYDA DTLKDPADG RKVKGVIIHWV SAAHALPVEI RLYDRLFSVP NPGAADDFLS VINPESLVIK QGFAEPSLKD AVAGKAFQFE REGYFCLDSR HSTAEPVFN RTVGLRDTWA KVGE |
| Source | E.coli |
| Target Names | glnS |
| Protein Names | Recommended name: Glutamine--tRNA ligase EC= 6.1.1.18 Alternative name(s): Glutaminyl-tRNA synthetase Short name= GlnRS |
| Expression Region | 2-554 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.