



# Recombinant Cysteine--tRNA ligase (cysS)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP813914EGX  |
| <b>Abbreviation</b>      | cysS   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q8FK44   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MLKIFNTLTR QKEEFKPIHA GEVGMVCGI TVYDLCHIGH GRTFVAFDVV<br>ARYLRFLGYK LKYVRNITDI DDKIIKRANE NGESFVALVD RMIAEMHKDF<br>DALNILRPDM EPRATHHIAE IIELTEQLIA KGHAYVADNG DVMFDVPTDP<br>TYGVLSRQDL DQLQAGARVD VVDDKRNPMDFVLWKMSKEG<br>EPSWSPWGA GRPGWHIECS AMNCKQLGNH FDIHGGGSDL<br>MFPHHENEIA QSTCAHDGQY VNYWMHSGMV MVDREKMSKS<br>LGNFFTVRDV LKYYDAETVR YFLMSGHYRS QLNYSEENLK QARAALERLY<br>TALRGTDKTV APAGGEAFAE RFIEAMDDDF NTPEAYSVLF DMAREVNRLK<br>AEDMAAANAM ASHLRKL SAV LGLLEQEPEA FLQSGAQADD SEVAEIEALI<br>QQRLDARKAK DWAAADAARD RLNEMGIVLE DGPQGTTWRR K |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | cysS   |
| <b>Protein Names</b>     | Recommended name: Cysteine--tRNA ligase EC= 6.1.1.16 Alternative name(s): CysteinyI-tRNA synthetase Short name= CysRS  |
| <b>Expression Region</b> | 1-461  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |