



Recombinant S-adenosylmethionine:tRNA ribosyltransferase-isomerase (queA)

Product Code	CSB-BP813916EGX
Abbreviation	queA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FKC5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MRVTD FSFEL PESLIAHYPM PERSSCRLLS LDGPTGALTH GTFTDLLDKL NPGDLLVFNN TRVIPARLFG RKASGGKIEV LVERMLDDKR ILAHIRASKA PKPGAELLG DDESINATMT ARHGALFEVE FNDQRSVLDI LNSIGHMPLP PYIDRPDEDA DRELYQTVYS EKPGAVAAPT AGLHFDEPLL EKLRAGVEM AFVTLHVGAG TFQPVVRVDI EDHIMHSEYA EVPQDVVDAV LAAKARGNRV IAVGTTSVRS LESAAQAAKN DLIEPLFDDT QIFIYPGFQY KVVDAVLTNF HLPESTLIML VSAFAGYQHT MNAYKAAVEE KYRFFSYGDA MFITYNPQAI NERVGE
Source	Baculovirus
Target Names	queA
Protein Names	Recommended name: S-adenosylmethionine:tRNA ribosyltransferase-isomerase EC= 5.-.-.- Alternative name(s): Queuosine biosynthesis protein QueA
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.