



Recombinant Protein fixA (fixA)

Product Code	CSB-BP813924EGX
Abbreviation	fixA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FLA3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKIITCYKCV PDEQDIANN ADGSLDFSKA DAKISQYDLN AIEAACQLKQ QAAEAQVTAL SVGGKALTNA KGRKDVLSRG PDELIVVIDD QFEQALPQQT ASVLAAAQK AGFDLILCGD GSSDLYAQQV GLLVGEILNI PAVNGVSKII SLTADTLTVE RELEDETETL SIPLPAVVAV STDINSPQIP SMKAILGAAK KPVQVWSAAD IGFNAEAAWS EQQVAAPKQR ERQRIVIEGD GEEQIAAF AE NLRKVI
Source	Baculovirus
Target Names	fixA
Protein Names	Recommended name: Protein fixA
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>