



# Recombinant *Oryza sativa* subsp. *japonica* Polygalacturonase inhibitor 1 (FOR1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP814080OFG   |
| <b>Abbreviation</b>      | FOR1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q8GT95  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | VRCPPS DKQALMRVKQ SLGNPATLST WSLASADCCE WDHVRCDEAG<br>RVNNVFIDGA NDVRGQIPSA VAGLTALMSL SLFRLPGLSG PIPACLTALS<br>NLQFLTISHT NVSGVIPDSL ARIRSLDSVD LSHNSLTGPI PNSFSDLPNL<br>RSLDLRSNKL TGCIPAGLVQ GQFRSLILSY NQLTGPIPRD DAQDEINTVD<br>LSHNRLTGDA SFLFAAGRPI GKVDLSWNDL DFDLSKLVFP PELTYLDLSH<br>NRIRGTVPRS LAALSTLQTL DLSYNRLCGP LPRLHGVIH GCKPYEHNQC<br>AGGAPLGGCH QS  |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | FOR1  |
| <b>Protein Names</b>     | Recommended name: Polygalacturonase inhibitor 1 Alternative name(s): Floral organ regulator 1 Polygalacturonase-inhibiting protein 1 Short name= PGIP-1   |
| <b>Expression Region</b> | 25-332  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |