



Recombinant Arabidopsis thaliana Tubby-like F-box protein 2 (TULP2)

Product Code	CSB-MP814086DOA
Abbreviation	TULP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GVE5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSLKSILRDL KEVRDGLGGI SKRSWSKSSH IAPDQTTPL DNIPQSPWAS LPPELLHDII WRVEESETAW PARAAVVSCA SVCKSWRGIT MEIVRIPEQC GKLTFPISLK QPGPRDSPIQ CFIKRN RATA TYILYYGLMP SETENDKLLL AARRIRRATC TDFIISLSAK NFSRSSSTYV GKLRSGFLGT KFTIYDNQTA SSTAQAQPNR RLHPKQAAPK LPTNSSTVGN ITYELNVLRT RGPRRMHCAM DSIPLSSVIA EPSVVQGIEE EVSSSPSPKG ETITTDKEIP DNPSLRDQP LVLKNKSPRW HEQLQCWCLN FKGRVTVASV KNFQLVAEID ASLDAPPEEH ERVILQFGKI GKDIFTMDYR YPLSAFQAF A ICISSEFDTKP ACEG
Source	Mammalian cell
Target Names	TULP2
Protein Names	Recommended name: Tubby-like F-box protein 2 Short name= AtTLP2
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.