



Recombinant Arabidopsis thaliana Probable isoaspartyl peptidase/L-asparaginase 2 (At3g16150)

Product Code	CSB-MP814094DOA
Abbreviation	At3g16150
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GXX1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGGWAIAVHG GAGIDPNLPA ERQEEAKQLL TRCLNLGIIA LRSNVSAIDV VELVIRELET DPLFNSSGRGS ALTEKGTVEM EASIMDGTKR RCGAVSGITT VKNPISLARL VMDKSPHSYL AFSGAEDFAR KQGVEIVDNE YFVTDDNVGM LKLAKANSI LFDYRIPPMG CAGAAATDSP IQMNGLPISI YAPE
Source	Mammalian cell
Target Names	At3g16150
Protein Names	Recommended name: Probable isoaspartyl peptidase/L-asparaginase 2 EC= 3.4.19.5 Alternative name(s): L-asparagine amidohydrolase 2 Cleaved into the following 2 chains: 1. Isoaspartyl peptidase/L-asparaginase 2 subunit alpha 2. Isoasp
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.