



# Recombinant *Macaca mulatta* C-C motif chemokine 5 (CCL5)

<b>Product Code</b>	CSB-YP814163MOW
<b>Abbreviation</b>	CCL5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8HYQ1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca mulatta</i> (Rhesus macaque)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SPHASDT TPCCFAYIAR PLPRAHIKEY FYTSGKCSNP AVVFVTRKNR QVCANPEKKW VREYINSLEM S
<b>Source</b>	Yeast
<b>Target Names</b>	CCL5
<b>Protein Names</b>	Recommended name: C-C motif chemokine 5 Alternative name(s): Small-inducible cytokine A5 T-cell-specific protein RANTES
<b>Expression Region</b>	24-91
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene is one of several CC cytokine genes clustered on the q-arm of chromosome 17. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene functions as a chemoattractant for blood monocytes, memory T helper cells and eosinophils. It causes the release of histamine from basophils and activates eosinophils. This cytokine is one of the major HIV-suppressive factors produced by CD8+ cells. It functions as one of the natural ligands for the chemokine receptor CCR5 and it suppresses in vitro replication of the R5 strains of HIV-1, which use CCR5 as a coreceptor.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.