



Recombinant Human Coiled-coil domain-containing protein 28A (CCDC28A)

Product Code	CSB-YP814216HU
Abbreviation	CCDC28A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IWP9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPRAEPRATL GEQEKAGLPL GAWRLYLLRH FRKQTELRRS GSRDVTGALL VAAAVASEAV GSLRVAEGGP NTLQLVLRW WPWCNKELKT MEERKVKRRS PKSFAHCTQ VVNAKKNAIP VSKSTGFSNP ASQSTSQRPK LKRVMKEKTK PQGGEGKGAQ STPIQHSFLT DVSDVQEMER GLLSLLNDFH SGKLQAFGNE CSIEQMEHVR GMQEKLARLN LELYGELEEL PEDKRKTASD SNLDRLLSDL EELNSSIQKL HLADAQDVPN TSAS
Source	Yeast
Target Names	CCDC28A
Protein Names	Recommended name: Coiled-coil domain-containing protein 28A Alternative name(s): CCRL1AP
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is located in a region close to the locus of the pseudogene of chemokine (C-C motif) receptor-like 1 on chromosome 6. The specific function of this gene has not yet been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.