



Recombinant Human Outer dense fiber protein 3-like protein 1 (ODF3L1)

Product Code	CSB-YP814225HU
Abbreviation	ODF3L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXM7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKLPKGTRSS VYFAQHPEKE PLPSRQEVKQ TPVIMAKIKG PGPAYLRLPS CTGYIDHDIS MFKAPAYTLH SRHSEKRMVC HSSPGPCYLL DPKITRFGMS SCPQVPMEEER ISNLRNPTL ASCQYYFEKI HPPGERRAPQ YTFGYRRPYR VMDLNPAPNQ YQMPLLLGPN TPVSRAAPCY SLASRDKNWF YKEDVAGGPG PTTYARPEPS IYQNRSPTY S MAKRFAYPLD LTPRPGPGSH EVQQVTVHKP HIPAFTMGIK HSLHLCPLVI DIRD
Source	Yeast
Target Names	ODF3L1
Protein Names	Recommended name: Outer dense fiber protein 3-like protein 1
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.