



Recombinant Human Protein FAM126B (FAM126B)

Product Code	CSB-MP814227HU
Abbreviation	FAM126B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXS8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLGTDRCVVE EWLSEFKALP DTQITSYAAT LHRKKTLPVA LYKVIQDSNN ELLEPVCHQL FELYRSSEVR LKRFTLQFLP ELMWVYLRLT VSRDRQSNCG IEALLLGIYN LEIADKDGNN KVLSTIPSL SKPSIYHEPS TIGSMALTEG ALCQHDLIRV VYSDLHPQRE TFTAQNRFEV LSFLMLCYNS AIVYMPASSY QSLCRMGSRV CVSGFPRQHE KHWKELCGRI VLDPEFMVQL LTGVYYAMYN GQWDLGQEV LDDIYRAQLE LFSQPLLVAN AMKNSLPFDA PDSTQEGQKV LKVEVTPTVP RISRTAITTA SIRRHRWRRE GAEGVNGGEE SVNLNDADEG FSSGASLSSQ PIGTKPSSSS QRGSLRKVAT GRSKDKETA SAIKSSSEPR DSVVRKQYVQ QPTDLSVDSV ELTPMKKHLS LPAGQVVPKI NSLSLIRTAS ASSKSFYV NGSQASTSIG VGTEGGTNLA ANNANRYSTV SLQEDRLGQA GEGKELLSPG APLTKQSRSP SFNMQLISQV
Source	Mammalian cell
Target Names	FAM126B
Protein Names	Recommended name: Protein FAM126B
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.