



Recombinant *Oryza sativa* subsp. *japonica* SAP-like protein BP-73 (BP-73)

Product Code	CSB-EP814457OFG-B
Abbreviation	BP-73
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8L4E7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MLLPPHPGRT LPAYHGDNRF LLGACFLSKL PMLRPMKLSL VCSANPNNHR SRSSDITRHQ KGGSARRKSK PYQEKDDSEN IDEFDTDIMS SKNGPPISLT SNSRPQATSV PGEREKEIVE LFKRVQAQLR ARGKGKEEKK PEQAKAQGER GSVDSLLNLL RKHSVDQRRK SGDEKEQSVD QTKRSNESGN KQNSSIFIKN DTQEEQKKPH PAAFKRPASN FRRRSPVNV KFPVTVNDA ERVINNINDA VQEAAPTLEN KAATDEPDSV STFEPNSVIE PENLSLDDLD HISDDEPDAS DTDEPSGEYD EPSLQIPVVP IIDESHTTL KSSLGGPDLs TLKVTELREL AKSRGIKGYs KMKKNDLVEL LSNMA
Source	E.coli
Target Names	BP-73
Protein Names	Recommended name: SAP-like protein BP-73 Alternative name(s): OsBP-73 Riaa1
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.