



# Recombinant Arabidopsis thaliana Protein BASIC PENTACYSTEINE6 (BPC6)

<b>Product Code</b>	CSB-MP814478DOA
<b>Abbreviation</b>	BPC6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8L999
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDDGGHRENG RHKAAVQGQW LMQHQPSMKQ VMSIIAERDA AIQERNLAIS EKKAAVAERD MAFLQRDTAI AERNNAIMER DSALTALQYR ENSMVTAPAA NMSACPPGCQ ISRGVKHLHH PHMHHHHQQH HIPQLTENAY ETREMEPNDG LPTSPAGST LESAKPKRGK RVNPKATTQT AANKRGPKNQ RKVKKESEDD LNKIMFVKTT HDYTDDEDSSK HILIGSKSDW KSQEMVGLNQ VVDYDETTMPP PVCSTGVLR QCYKWGNGGW QSSCCTTTL MYPLPALPNK RHARVGGGRKM SGSAFNKLLS RLAAEGHHDL SNPVDLKDHW AKHGTRNYIT IK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BPC6
<b>Protein Names</b>	Recommended name: Protein BASIC PENTACYSTEINE6 Short name= AtBPC6 Alternative name(s): GAGA-binding transcriptional activator BBR/BPC6
<b>Expression Region</b>	1-342
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.