



Recombinant Arabidopsis thaliana Methyl-CpG-binding domain-containing protein 2 (MBD2)

Product Code	CSB-EP814482DOA
Abbreviation	MBD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LA53
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSMSQSRAVQ RSSSPNEDRG ENQLVVYDLK GNDDTEEEVL PVQSQPLSSR TQCPSIGAFT VQCASCFKWR LMPSMQKYEE IREQLLNPF FCDTAREWKP DISCDVPADI YQDGTRLWAI DKPNISRPPA GWQRLLRIRG EGGTRFADV YVAPSGKKLR STVEVQKYLN DNSEYIGEGV KLSQFSFQIP KPLQDDYVRK RPARLLDSID NTNTPVAKEA NPLAWISPDH HSLQLGTPT ESGLNNSHYQ PSKKKKTSTL SIFGSNDELA DR
Source	E.coli
Target Names	MBD2
Protein Names	Recommended name: Methyl-CpG-binding domain-containing protein 2 Short name= AtMBD2 Short name= MBD02 Alternative name(s): Methyl-CpG-binding protein MBD2
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.