



# Recombinant Mouse CAP-Gly domain-containing linker protein 4 (Clip4)

<b>Product Code</b>	CSB-EP814551MO-B
<b>Abbreviation</b>	Clip4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8CI96
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTIEDLPDIP LEGSSLIGRY PFLFTGSDTS VIFSISAAPM PSDCEFSFFD PNDASCQEIL FDPKTSVSEL FAILRQWVPQ VQQNIDIIGN EILKRGCVN DRDGLTDMTL LHYTCKSGAH GIGDIETAVK FAAQLIDLGA DASLRSRWTN MNALHYASYF DVPELIRVIL KTSKPKDVDA TCSDFNFGTA LHIAAHNLCA GAVKTLLELG ANPAFRNDKG QIPAEVVPDP VDMPLEMADA AAIAKEIKQM LLDAVPLPCT SAKAVLPNSD HTTSRAMLTS LGLKLGDRVV IAGQKVGTLR FCGTTEFASG QWAGIELDEP EGKNNGSVGR VQYFKCAPKY GIFAPLSKIS KLDGGRKTTT HTPSTRATPH ARSQKVDVAH FTSRVNSGLT TSKKETASES TLTLPPSEEP KTVAENDAAQ PGSMSSSSSS SSLDHKQSYK KKLTTSSGGK KTLKSPSLP SRASAGLKSS ATSAANNSHH EGALHLGERV LVVGQRVGTI KFFGTTNFAP GYWYGIELEK PHGKNDGSVG GVQYFSCSPR YGIFAPPSRV QRLSDSLDTL SEISSNKQNH SYPGFRRSFS TTSASSQKEI NRRNAFAKTK TTLRRSWSSS TTAGGLEGTV KLHEGSQVLL TSSNEMATVR YVGPTDFASG IWLGLELRSA KGKNDGAVGD KRYFTCKPNY GVLVPRSRVT YRGISGSKLI DENC
<b>Source</b>	E.coli
<b>Target Names</b>	Clip4
<b>Protein Names</b>	Recommended name: CAP-Gly domain-containing linker protein 4 Alternative name(s): Restin-like protein 2
<b>Expression Region</b>	1-704
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .