



Recombinant sn-glycerol-3-phosphate-binding periplasmic protein ugpB (ugpB)

Product Code	CSB-BP814662EGX
Abbreviation	ugpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8CVL9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	VTTIPFW HSMEGELGKE VDSLARFNA ENPDYKIVPT YKGNIEQNLS AGIAAFRTGN APAILQVYEV GTATMMASKA IKPVYDVFKE AGIQFDESQF VPTVSGYYSD SKTGHLLSQP FNSSTPVLYY NKDAFKKAGL DPEQPPKTWQ DLADYAAKLK ASGMKCGYAS GWQGWIQLEN FSAWNGLPFA SKNNGFDGTD AVLEFNKPEQ VKHIAMLEEM NKKGDFSYVG RKDESTTEKPY NGDCAMTTAS SGSLANIREY AKFNYGVGMM PYDADAKDAP QNAIIGGASL WVMQGKDKET YTGVAKFLDF LAKPENAAEW HQKTGYLPIT KAAIDLTTREQ GFYEKNPGAD IATRQMLNKP PLPFTKGLRL GNMPQIRVIV DEELESVWTG KKTPQQALDT AVERGNQLLR RFEKSTKS
Source	Baculovirus
Target Names	ugpB
Protein Names	Recommended name: sn-glycerol-3-phosphate-binding periplasmic protein ugpB
Expression Region	24-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.