



Recombinant Phosphomethylpyrimidine synthase (thiC)

Product Code	CSB-MP815184EGX
Abbreviation	thiC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FB77
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSATKLTRRE QRAQAQHFID TLEGSAFPNS KRIYITATQP GVRVPMREIQ LSPTLIGGSK EQPQYEENE IPVYDTSGPY GDPQIAINVQ QGLAKLRQPW IDARGDTEEL TVRSSDYTKA RLADDGLDEL RFGVLT PKR AKAGRRVTQL HYARQGIITP EMEFIAIREN MGRERIRSEV LRHQHPGMSF GARLPENITA EFVRDEVAAG RAIIPANINH PESEPMIIGR NFLVKVNANI GNSAVTSSIE EEVEKLVWST RWGADTVMDL STGRYIHETR EWILRNVPV IGTVPYIYQAL EKVNGIAEDL TWEVFRDILL EQAEQGVDFY TIHAGVLLRY VPMTAKRLTG IVSRGGSIMA KWCLSHHQEN FLYQHFREIC EICAAYDVSL SLGDGLRPGS IQDANDEAQF AELHTLGELT KIAWEYDVQV MIEGPGHVPM QMIRRNMTTE LEHCHEAPFY TLGPLTTDIA PGYDHFTSGI GAAMIGWFGC AMLCYVTPKE HLGLPNKEDV KQGLITYKIA AHAADLAKGH PGAQIRDNAM SKARFEFRWE DQFNALDPF TARAYHDETL PQESGKVAHF CSMCGPKFCS MKISQEV RDY AAAQTIEVGM ADMSENFRRAR GGEIYLRKEE A
Source	Mammalian cell
Target Names	thiC
Protein Names	Recommended name: Phosphomethylpyrimidine synthase EC= 4.1.99.17 Alternative name(s): Hydroxymethylpyrimidine phosphate synthase Short name= HMP-P synthase Short name= HMP-phosphate synthase Short name= HMPP synthase
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.