



Recombinant Probable UDP-N-acetyl-D-mannosaminuronic acid transferase (wecG)

Product Code	CSB-MP815189EGX
Abbreviation	wecG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FBP7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNNNTTAPTY TLRGLQLIGW RDMQHALDYL FADGQLKQGT LVAINAEKML TIEDNAEVRE LINAAEFKYA DGISVVRSVR KKYPQAQVSR VAGADLWEEL MARAGKEGTP VFLVGGKPEV LAQTEAKLRN QWNVNIVGSQ DGYFKPEQRQ ALFERIHASG AQIVTVAMGS PKQEIFMRDC RLVHPDALYM GVGTTYDVFT GHVKRAPKIW QTLGLEWLYR LLSQPSRIKR QLRLRLRYLRW HYTGNL
Source	Mammalian cell
Target Names	wecG
Protein Names	Recommended name: Probable UDP-N-acetyl-D-mannosaminuronic acid transferase Short name= UDP-ManNAcA transferase EC= 2.4.1.-
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.