



Recombinant L-seryl-tRNA (Sec) selenium transferase (selA)

Product Code	CSB-MP815194EGX
Abbreviation	selA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FCC1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MTTETRSLYS QLPAIDRLLR DSSFLSLRDT YGHTRVVELL RQMLDEAREV IRDSQTLPAW CENWAQEVDA RLTKEAQSAL RPVINLTGTV LHTNLGRALQ AEAAVEAVTK AMRSPVTLEY DLDDAGRGRH DRALAQLLCR ITGAEDACIV NNNAAVLLM LAATASGKEV VVSRGELVEI GGAFRIPDVM RQAGCTLHEV GTTNRTHAND YRQAVNENTA LLMKVHTSNY SIQGFKAID EAELVALGKE LDVPVVDLG SGSLVDLSQY GLPKEPMPQE LIAAGVSLVS FSGDKLLGGP QAGIIVGKKE MIARLQSHPL KRALRADKMT LAALATLRL YLHPEALSEK LPTLRLLTRS AEVIQIQAGR LQAPLAHYG AEFVQVMPC LSQIGSGSLP VDRLPSAALT FPHDGRGSH LESLAARWRE LPVPVIGRIY DGRLWDLRC LEDEQRFLEM LLK
Source	Mammalian cell
Target Names	selA
Protein Names	Recommended name: L-seryl-tRNA(Sec) selenium transferase EC= 2.9.1.1 Alternative name(s): Selenocysteine synthase Short name= Sec synthase Selenocysteinyl-tRNA(Sec) synthase
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.