



# Recombinant Glutamate--tRNA ligase (gltX)

<b>Product Code</b>	CSB-BP815218EGX
<b>Abbreviation</b>	gltX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FFC9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIKTRFAPS PTGYLHVGG A RTALYSWLFA RNHGGEFVLR IEDTDLERST PEAIEAIMDG MNWLSLEWDE GPYYQTKRFD RYNAVIDQML EEGTAYKCYC SKERLEALRE EQMAKGEKPR YDGRCRHSHE HHADDEPCVV RFANPQEGSV VFDDQIRGPI EFSNQELDDL IIRRTDGSPT YNFCVVVDDW DMEITHVIRG EDHINNTPRQ INILKALKAP VPVYAHVSMI NGDDGKKLSK RHGAVSVMQY RDDGYLPEAL LNYLVRLGWS HGDQEIFTRE EMIKYFTLNA VSKSASAFNT DKLLWLNHHY INALPPEYVA THLQWHIEQE NIDTRNGPQL ADLVKLLGER CKTLKEMAQS CRYFYEDFAE FDADAACKHL RPVARQPLEV VRDKLTAITD WTAENVHAI QATADELEVG MGKVGMPLRV AVTGAGQSPA LDVTVHAIGK TRSIERINKA LAFIAERENQ Q
<b>Source</b>	Baculovirus
<b>Target Names</b>	gltX
<b>Protein Names</b>	Recommended name: Glutamate--tRNA ligase EC= 6.1.1.17 Alternative name(s): Glutamyl-tRNA synthetase Short name= GluRS
<b>Expression Region</b>	1-471
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.