



Recombinant Hydroxyethylthiazole kinase (thiM)

Product Code	CSB-EP815224EGX-B
Abbreviation	thiM
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FFY0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQVDLLSSAQ SAHALHLFHQ HSPLVHCMTN DVVQTFNTANT LLALGASPAM VIETEEASQF AAIASALLIN VGTLTQPRAQ AMSAAVEQAT RSQTPWTLDP VAVGALDYRR RFCVELLSHK PTAIRGNASE IMALAGVANG GRGVDTTDAA ANAIPAAQTL ARETGAIVVV TGEVDYVTDG HRIIGIHGGD PLMTKVVGTG CALSAVVAAC CALPGDTLEN IASACHWMKQ AGERAVARSE GPGSFVPHFL DALWQLAQEV QA</p>
Source	E.coli
Target Names	thiM
Protein Names	<p>Recommended name: Hydroxyethylthiazole kinase EC= 2.7.1.50 Alternative name(s): 4-methyl-5-beta-hydroxyethylthiazole kinase Short name= TH kinase Short name= Thz kinase</p>
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>