



Recombinant Hydroxyethylthiazole kinase (thiM)

Product Code	CSB-YP815224EGX
Abbreviation	thiM
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8FFY0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQVDLLSSAQ SAHALHLFHQ HSPLVHCMTN DVVQTFTANT LLALGASPAM VIETEEASQF AAIASALLIN VGTLTQPRAQ AMSAAVEQAT RSQTPWTLDP VAVGALDYRR RFCVELLSHK PTAIRGNASE IMALAGVANG GRGVDTTDAA ANAIPAAQTL ARETGAIVVV TGEVDYVTDG HRIIGIHGGD PLMTKVVG TG CALSAVVAAC CALPGDTLEN IASACHWMKQ AGERAVARSE GPGSFVPHFL DALWQLAQEV QA</p>
Source	Yeast
Target Names	thiM
Protein Names	<p>Recommended name: Hydroxyethylthiazole kinase EC= 2.7.1.50 Alternative name(s): 4-methyl-5-beta-hydroxyethylthiazole kinase Short name= TH kinase Short name= Thz kinase</p>
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>