



# Recombinant Protein ydeP (ydeP)

<b>Product Code</b>	CSB-EP815235EGX
<b>Abbreviation</b>	ydeP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8FHF8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKKKIESYQG AAGGWGAVKS VANAVRKQMD IRQDVAMFD MNKPEGFDCP GCAWPDPKHS ASFDICENGA KAIWEVTDK QVNASFFAEN TVQSLLTWGD HELEAAGRLT QPLKYDDVSD CYKPLSWQQA FDEIGARLQS YSDPNQVEFY TSGRTSNEAA FLYQLFAREY RSNNFPDCSN MCHEPTSVGL AASIGVGKGT VLLEDFEKCD LVICIGHNPG TNHPRMLTSL RALVKRGAKM IAINPLQERG LERFTAPQNP FEMLTNSETQ LASAYYNVRI GGDMALLKGM MRLLIERDDA ASAAGRPSLL DDEFIQHTTV GFDELRRDVL NSEWKDIERI SGLSQTQIAE LADAYAAAER TIICYGMGIT QHEHGTQNVQ QLVNLLLMMKG NIGKPGAGIC PLRGHSNVQG DRTVGITEKP SAEFLDRLCE RYGFTPPHAP GHAAIASMQA ICTGQARALI CMGGNFALAM PDREASAVPL TQLDLAVHVA TKLNRSHLLT ARHSYILPVL GRSEIDMQKS GAQAVTVEDS MSMIHASRGV LKPAGVMLKS ECAVVAGIAQ AALPQSVVAW EYLVEDYDRI RNDIEAVLPE FADYNQRIRH PGGFHLINAA AERRWMTSSG KANFITSKGL LEDPSSAFNS KLVMATVRSH DQYNTTIYGM DDRYRGVFGQ RDVVFMSAKQ AKICRVKNGE RVNLIALTPD GKRSSRRMDR LKVVIYPMAD RSLVTYFPES NHMLTLDNHD PLSGIPGYKS IPVELEPSN
<b>Source</b>	E.coli
<b>Target Names</b>	ydeP
<b>Protein Names</b>	Recommended name: Protein ydeP
<b>Expression Region</b>	1-759
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.