



Recombinant Probable malonic semialdehyde reductase RutE (rutE)

Product Code	CSB-EP815244EGX-B
Abbreviation	rutE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FJ44
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNEAVSPGAL STLFTDARTH NGWRETPVSD ETLREIYALM KWGPTSANCS PTRIVFIRTA EGKERLRPAL SSGNLQKTLT APVTAIVAWD SEFYERLPQL FPHGDARSWF TSSPQLAEET AFRNSSMQAA YLIVACRALG LDTGPMSGFD RQYVDDAFFA GSTLKSNLLI NIGYGDSSKL FARLPRLSFE EACGLL
Source	E.coli
Target Names	rutE
Protein Names	Recommended name: Probable malonic semialdehyde reductase RutE EC= 1.1.1.298
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.