



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At3g60040 (At3g60040)

Product Code	CSB-EP815427DOA-B
Abbreviation	At3g60040
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GW12
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDAKMFDTIC SRDAISWLPD EVLGKILSLI PTKQAVSTSL LAKKWRTIFR LVDHLELDDDS FSLQAVKDQT PGLRKHVRFV FTEDFKIFVD RTLALQCDYP IKNFSLKCHV SKYDERQKAC VGRWISNVVG RGVFELDLRM KDPGIHFLPP HLVASKTLVK LTLGTQLCLG QLPSYVSLPS LKSLFIDTIV FYDIEDLCCV LLAGCPVLEE LSVHHDFIA TPHTISSPTL KRLSVDYHCP DDVDSASHMS FDLPKLVYLE YSHCALGEYW QINLESLVEA KLDLGLERKV LRMDVTDLII GIRNVQSLHL SPDSVHVIYS YCRHGLPVFK NLVNLSFGSK NKRGRWLLAN LLKQSTKLET LIVKDLNGYT GDVSMPLNKV TSLHILGYRG TADEVKQLKS FIGEFECLEL VQVDVAEAAE DNGKILQSKR DLMMLLGVSL PSKCQFKVT
Source	E.coli
Target Names	At3g60040
Protein Names	Recommended name: F-box/LRR-repeat protein At3g60040
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.