



Recombinant Arabidopsis thaliana 1,2-dihydroxy-3-keto-5-methylthiopentene dioxygenase 4 (ARD4)

Product Code	CSB-YP815453DOA
Abbreviation	ARD4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H185
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MALEAWFMDD SNEDQRLPHH RNPKELVSLD YLAELGVLYW KLNPENYEND SELSKIREDR GYDYMDLLDL CPEKVSNYEE KLKNFFTEHI HKDEEIRYCL AGSGYFDVRD KDDRWIRIWM QPGDLIVLPA GIYHRFTLDA SNYIKLMRLF VGEPVWTPYN RPQEEHPVRK KYIHGLTYKF GETVKAH
Source	Yeast
Target Names	ARD4
Protein Names	Recommended name: 1,2-dihydroxy-3-keto-5-methylthiopentene dioxygenase 4 EC= 1.13.11.54 Alternative name(s): Acireductone dioxygenase (Fe(2+)- requiring) 4 Short name= ARD 4 Short name= Fe-ARD 4
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.