



Recombinant Arabidopsis thaliana 1-aminocyclopropane-1-carboxylate oxidase homolog 3 (At1g06650)

Product Code	CSB-MP815456DOA
Abbreviation	At1g06650
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8H1S4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEMMKIDPLF DRASELKAFD ETKTGVKGLV DSGVSVQVPRI FHHPTVKLST PKPLPSDLLH LKTIPTIDLG GRDFQDAIKR NNAIEEIKEA AAKWGGFFQVI NHGVSLELLE KMKKGVRDFH EQSQEVRKEF YSRDFSRRFL YLSNFDLFFSS PAANWRDTFS CTMAPDTPKP QDLPEICRDI MMEYSKQVMN LGKFLFELLS EALGLEPNHL NDMDCSKGLL MLSHYPPCP EPDLTLGTSQ HSDNSFLTIVL LPDQIEGLQV RREGHWFDVP HVSGALIINI GDLLQLITND KFISLEHRVL ANRATRARVS VACFFTTGVR PNPRMYGPIR ELVSEENPPK YRETTIKDYA TYFNAKGLDG TSALLHFKI
Source	Mammalian cell
Target Names	At1g06650
Protein Names	Recommended name: 1-aminocyclopropane-1-carboxylate oxidase homolog 3
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.