



# Recombinant Human Interferon regulatory factor 2-binding protein 1 (IRF2BP1)

<b>Product Code</b>	CSB-BP815541HU
<b>Abbreviation</b>	IRF2BP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IU81
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASVQASRRQ WCYLCDLPKM PWAMVWDFSE AVCRCGVNFE GADRIELLID AARQLKRSHV LPEGRSPGPP ALKHPATKDL AAAAAQGPQL PPPQAQPQPS GTGGGVSGQD RYDRATSSGR LPLPSALEY TLGSRLANGL GREEAVAEGA RRALLGSMPG LMPPGLLAAA VSGLGSRGLT LAPGLSPARP LFGSDFEKEK QQRNADCLAE LNEAMRGRAE EWHGRPKAVR EQLLALSACA PFNVRFKKDH GLVGRVFAFD ATARPPGYEF ELKLFTEYPC GSGNVYAGVL AVARQMFHDA LREPGKALAS SGFKYLEYER RHGSGEWRQL GELLTDGVRS FREPAPAEAL PQQYPEPAPA ALCGPPPRAP SRNLAPTPRR RKASPEPEGE AAGKMTTEEQ QQRHWVAPGG PYSAETPGVP SPIAALKNVA EALGHSPKDP GGGGGPVRAG GASPAASSTA QPPTQHRLVA RNGEAEVSPT AGAEAVSGGG SGTGATPGAP LCCTLCRERL EDTHFVQCPS VPGHKFCFPC SREFIKAQGP AGEVYCPSGD KCPLVGSSVP WAFMQGEIAT ILAGDIKVKK ERDP
<b>Source</b>	Baculovirus
<b>Target Names</b>	IRF2BP1
<b>Protein Names</b>	Recommended name: Interferon regulatory factor 2-binding protein 1 Short name= IRF-2-binding protein 1 Short name= IRF-2BP1
<b>Expression Region</b>	1-584
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.