



Recombinant Human Keratin-associated protein 13-1 (KRTAP13-1)

Product Code	CSB-MP815543HU
Abbreviation	KRTAP13-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IUC0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYNCCSGNF SSRSCGGYLH YPASSCGFSY PSNQVYSTDL CSPSTCQLGS SLYRGCQQTC WEPTSCQTSY VESSPCQTSC YRPRTSLLCS PCQTTYSGSL GFGSSSCRSL GYGSRSYSV GCGSSGFRSL GYGCGFPSL GYGVGFCRPT YLASRSCQSS CYRPTCGSGF YY
Source	Mammalian cell
Target Names	KRTAP13-1
Protein Names	Recommended name: Keratin-associated protein 13-1 Alternative name(s): High sulfur keratin-associated protein 13.1
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.