



# Recombinant Human CREB3 regulatory factor (CREBRF)

<b>Product Code</b>	CSB-BP815547HU
<b>Abbreviation</b>	CREBRF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IUR6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPQPSVSGMD PPFGD AFRSH TFSEQ TLMST DLLANSSDPD FMYELDREM N YQQNPRDNFL SLEDCKDIEN LESFTDVL DN EGALTSNWEQ WDTYCEDLTK YTKLTSCDIW GTKEVDYLGL DDFSSPYQDE EVISK TPTLA QLNSEDSQSV SDSLYYPDSL FSVKQNPLPS SFPGKKITSR AAAPVCSSKT LQAEVPLSDC VQKASKPTSS TQIMVKTNMY HNEKVN FHV E CKDYVKKAKV KINPVQQSRP LLSQIH TDAA KENTCYCGAV AKRQEKKGME PLQGHATPAL PFKETQELL SPLPQEGPGS LAAGESSLS ASTSVSDSSQ KKEEHNYSLF VSDNLGEQPT KCSPEEDEED EEDVDDDEDHD EGFGESEHLS ENEEEEEEEEE DYEDDKDDDI SDTFSEPGYE NDSVEDLKEV TSISSRKR GK RRYFWEYSEQ LTPSQQERML RPSEWNRDTL PSNMYQKNGL HHGKYAVKKS RRTDVEDLTP NPKKLLQIGN ELRKLNKVIS DLTPVSELPL TARPRSRKEK NKLASRACRL KKKAQYEANK VKLWGLNTEY DNLLFVINSI KQEIVNRVQN PRDERGPNMG QKLEILIKDT LGLPVAGQTS EFNQVLEKT AEGNPTGGLV GLRIPTSKV
<b>Source</b>	Baculovirus
<b>Target Names</b>	CREBRF
<b>Protein Names</b>	Recommended name: CREB3 regulatory factor Alternative name(s): Luman recruitment factor Short name= LRF
<b>Expression Region</b>	1-639
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.